



PRODUCT SPECIFICATION



Product: CEDERROTH EYE WASH

Reference code: 7251, 725200



1. General description

A personal eyewash bottle containing sterile, buffered, saline eye wash solution.

The bottle is easily opened and ready to use by twisting off the dustcover. The integrated eye cup helps keeping the eye open during rinsing and ensures a generous flow without unnecessary waste

2. Intended use

The Cederroth Eye Wash is intended for washing away foreign particles and dangerous liquids, like alkali or acids, from an eye in an emergency first aid situation at work places. The Cederroth buffered Eye Wash also has a neutralising effect on splashes from alkalis and acids. The effect is in general stronger for alkalis than for acids.

The Cederroth Eye Wash is intended for the first acute phase at the very accident site before getting to a plumbed or self-contained eye wash equipment or seeking medical attention. The Cederroth Eye Wash can also be used during transport to medical care.

The durable consumer box contains 2 Cederroth Eye Wash bottles.

3. Technical description

Product dimensions :

Bottle

Diameter (cm): 6.6
Height (cm) 23,5
Fill volume (ml): 500

Box (2 bottles)

l x w x h (cm): 7 x 14,5 x 24,5

Material

Liquid

Sterile, isotonic Sodium Chloride 0.9%,
Buffered with Sodium Borate to neutral PH.
Preservative free.

Bottle

Leak proof polypropylene bottle produced with Blow-Fill-Seal technique.
Eye cup in transparent polypropylene
Dust cover in green hard plastics (ABS)

4. Sterilization

Sterilization method: Aseptic filling



PRODUCT SPECIFICATION



5. Composition

Consumer package (725200)

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
7251	Cederroth Eye Wash	2

Transport package

Number of consumer packages in transport package:	6
Dimensions (l x w x h) mm:	300 x 235 x 260
Weight (kg):	7,8

6. Labeling

On product and consumer package

Cederroth address, telephone- and fax-number, e-mail and web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code
Instructions for use
Symbols: CE +NB#, STERILE A, DO NOT REUSE, LOT, USE BY, TEMPERATURE LIMITATION

On transport package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code
Quantity
Symbols: CE +NB#, LOT, USE BY, FRAGILE

7. Storage conditions

Store in temperatures between +2°C and +40°C.

Shelf life: 4,5 years from production date

8. Accessories

The following accessories can be used with the product but are not necessary.

7200

Wall bracket for Cederroth Eye Wash



For organized storage and defined placement. Also works as an anti-theft system

721500

Cederroth Eye Wash station



For storage of 2 Cederroth Eye Wash bottles. Plaster dispenser and empty compartment for additional products.

790400

Heating Cabinet for Cederroth Eye Wash



Authorization: