



PRODUCT SPECIFICATION



Product: CEDERROTH BURN GEL
DRESSING
Reference code: 901900



1. General description

A foam dressing impregnated with burn gel for the immediate first aid treatment of minor burns. The gel rapidly cools and soothes the burn injury, prevents further tissue damage and offers a rapid and ongoing pain relief.

The dressing is sterile and individually packed in a sealed pouch which easily can be teared open. Two Cederroth Burn Gel Dressings are packed in one consumer box.

The Cederroth Burn Gel Dressing fits into Cederroth First Aid Stations / First Aid Panels.

2. Intended use

The intended use of Cederroth Burn Gel Dressing is first aid treatment of 1st and 2nd degree of minor burns, mosquito bites and sunburns.

3. Technical description

Product dimensions:

<u>Pouch</u>		<u>Box (incl. 2 dressings)</u>	
Height (mm):	124	l x w x h (mm):	125 x 35 x 65
Width (mm)	134		
Dressing (mm)	100 x 100		

Ingredients/ Material

Gel	Melaleuca Oil 1% (Tea Tree Oil), > 90% water based gel
Dressing	An open cell polyester foam, 3 mm thick
Pouch	Laminated plastic/Al-foil/plastic
Box	Carton

4. Sterilization

Sterilization method: Gamma irradiation, minimum dosage 25 kGy

5. Composition

Consumer package

<u>Ref#</u>	<u>Product/Material</u>
901900	Cederroth Burn Gel Dressing

Pcs
2

Transport package

Number of consumer packages
in transport package:

24

Dimensions (l x w x h) mm:

260 x 205 x 180

Weight (kg):

24



PRODUCT SPECIFICATION



6. Labeling

On product and consumer package

Manufacturing address, Cederroth address, telephone- and fax-number, e-mail and web-address

Product name, Reference code

Lot-number, Expiry date, EAN-code

Instructions for use including warning text.

Symbols: CE NB#, STERILE, R, DO NOT USE IF PACKAGE IS DAMAGED (only on product not on consumer package), DO NOT REUSE, LOT, USE BY, TEMPERATURE LIMITATION

On transport package

Manufacturing address, Cederroth address, telephone- and fax-number, web-address

Product name, Reference code

Lot-number, Expiry date, EAN-code, Quantity.

Symbols: CE NB#, LOT, USE BY

7. Storage conditions

Store in temperatures between +5°C (41°F) and +30°C (86°F).

Shelf life: 5 years from production date

8. Accessories

No accessories are necessary for the use of this product

Authorization: 20110802/ QA&RA Mgr Wound Care, Cederroth AB